





www.uniraymedical.com

# THE UNIRAY PROMISE

## PROTECTION FOR THE PEOPLE WHO PROTECT LIFE

With over 3 decades of global experience in x-ray and radiation protection, we understand your role as a medical professional along with the risks you take to fulfill it. We know that you deserve the greatest level of protection while enabling efficiency, comfort, and style.

UniRay is your destination for radiation protection products that represent a synergy of technological innovation, applied research and dedicated care. Our products meet and exceed the highest safety standards of the industry.

We design protective apparel used in interventional radiology and cardiology, Cath labs, x-ray assisted surgery and general radiology. UniRay offers a complete range of products that are designed with both function and style in mind. Quality, weight, ergonomics and aesthetics form the core attributes of every UniRay product.





## UNIRAY EQUALS SAFETY

When you use UniRay products, you are guaranteed optimum protection against scattered radiation.

UniRay products are tested according to IEC 61331-1:2014/ DIN EN 61331-1:2016, which establishes the protection of aprons across a specified range of tube voltages from 60 kV to 150kV. The new standard IEC EN 61331-1:2014/DIN EN 61331-1:2016, adopted a new testing methodology called "Broad Beam Geometry" (BBG) method to measure all Radiation Protection Aprons including lead, lead-free and lead composite protection materials. This testing methodology determines the lead equivalence of protection material considering fluorescent radiation as well as backscatter by the material itself.

# **THE X-FACTOR**

Our radiation protective sheetings are manufactured for perfect uniformity, suppleness and durability.

# **LITELEAD**<sup>™</sup>

**Max**Lite<sup>™</sup>

LITELead<sup>™</sup> is our classic leaded material with proprietary advance binding system. LEV 0.25mm Pb - weight 3.40 kg/sqm

**KryptoLite**<sup>™</sup> with its multi-layer technology eliminates fluorescence radiation and back scatter thus providing optimum safety to the user. **KryptoLite**<sup>™</sup> is at least 15-20% lighter than any other lead-free material in the world. This environmentally friendly lead-free material is formulated with a mixture of bismuth and antimony. LEV 0.25mm Pb - weight 2.80 kg/sqm

**LiteGreen™** This Environmentally friendly lead-free material is formulated with a mixture of Bismuth, Tungsten and Antimony. LiteGreen<sup>™</sup> is the material of choice for users. LEV 0.25mm Pb - weight 3.10 kg/sqm

# **KRYPTOLITE**

LiteGreen™

**MaxLite**<sup>™</sup> uses the substitution technique over a large part of lead in protective sheeting to result in a product that is 20-30% lighter than most composite materials. LEV 0.25mm Pb - weight 3.10 kg/sqm



**Optixx**<sup>™</sup>

• Significant overlap on the front opening of

the vest as well as the

area) also provides

added protection

• Optimised weight

• Full overlap on skirt

skirt (covering the waist

distribution to reduce load

on back and shoulders

circumference for added

lower body protection

## Maxx<sup>™</sup>

- Complete overlap on the front opening of the vest, as well as the overlap of the skirt for added protection
- Optimised weight distribution to reduce load on the back and shoulders
- Full overlap on skirt circumference for added lower body protection



\09

# **02** DOUBLE-SIDED APRON

ALL ROUND PROTECTION SOLUTION



Aprons are available in S/M/L/XL and custom sizes Aprons are available in LEV 0.25, 0.35, 0.5 mm Pb

- Significant overlap in front for added protection
- Velcro closure for a perfect fit
- Flexible belt used as standard option, to release load on back and shoulder

## Maxx<sup>™</sup>

- Full overlap in front for total protection
- Flexible belt used as standard option, to release load on back and shoulder
- Smart side closure to prevent accidental exposure





\11

# **03** BACKSAVER APRON

### **UPTO 50% REDUCTION IN LOAD ON BACK**

UniRay's **Back**SAVER lifts the entire back panel of the apron with use of a membrane. An upthrust is created when pressure is exerted at specific points and directions. The thrust opposes the weight of the apron and sufficiently lifts it, taking the load off the shoulder, back and spine.

Available in Skirt & Vest, Double-Sided and Coat Apron Also available in magnetic closure.





• Floatation on shoulder







11

• Built in belt, easy to clean and disinfect









Complete frontal protection

- Smart and flexible insert with velcro closure
- Built-in shoulder relief pads





## Complete frontal protection

- Easy drop off from under a surgeon's scrubs when removing
- Optimised weight distribution



\15

## **06** ZIPPER APRON **VEST WITH ZIP CLOSURE**




## **07** MEASUREMENTS **PERSONAL X-RAY PROTECTION - HOW TO MEASURE**

These are body measurements; apparel size can be determined based on these measurements.


The length of the X-ray protection clothing should be selected in accordance with the International Standard IEC 61331-3:2014 with coverage to below the knee for complete protection.

### **ALL DIMENSIONS IN CM**



• Easy to clean for better hygiene





LENGTH OF APRON

**VEST LENGTH** 

CHEST CIRCUMFERENCE

WAIST CIRCUMFERENCE

HIP CIRCUMFERENCE

SKIRT LENGTH

# **08** THYROID SHIELDS & COLLARS

### THE SENSITIVE GLAND PROTECTION

All thyroid shields can be affixed to the apron using a push button.



#### BALANCE

Unique design to protect the thyroid & sternum



#### UNIVERSAL

An all time favourite

• Magnetic closure




• Velcro closure





## ELEGANCE

For a snug fit enabling total protection





## All round thyroid protection



\19

# **09** LOWER BODY PROTECTION

PATIENT APRON





WITH SNAP BUCKLE & VELCRO CLOSURE

LEV 0.5 mm Pb in set of 4 sizes



WITH METALLIC RING CLOSURE LEV 0.5 mm Pb in a set of 4 sizes





NETO FLEXIBLE GONAD SHIELD

LEV 1.0 mm Pb in set of 4 sizes



**GONAD/OVARIAN SHIELD** LEV 1.0 mm Pb in set of 3 sizes

> Available in universal size Available in 0.25, 0.35, 0.5 mm Pb


### **FEZ** Protection of frontal cortex

Top view



## **11** CUSTOMISATION **EXPRESS YOUR UNIQUE WAY**

We maintain a vast design bank which can be embroidered on pockets or printed by our in-house precision embroidery equipment. We also provide a name tag of doctors or of a hospital for easy attachment on the apron pocket.



## **12** SHOULDER PROTECTION **SLEEVES**

Shoulder sleeves significantly reduce radiation on the shoulder joint. Radiation through armholes is also reduced. UniRay's shoulder sleeves allow complete freedom of hand movement during procedures. Sleeves are available in both Sewn-on and Button-on versions in either only left, only right or both shoulders versions.





### **SEWN-ON SLEEVE**

- Factory attached
- Also available in full sleeves

#### **BUTTON-ON SLEEVE**


- Attaches easily to all models of Vests and Double-Sided Aprons
- Also available in full sleeves

## **13** FABRICS AND COLORS VIBRANT COLORS THAT LIFT MOODS


VIDRANT COLORS THAT LIFT WOODS

 $\mathsf{UNITEX^{\mathsf{TM}}},$  is antibacterial, comfortable & stain resistant.


It is Invented as a soft & anti-wrinkle covering for all radiation protection apparel.




**FUSCHIA** Joyful



ROYAL BLUE Energetic



**LEMON** Cheerful




**MOUNTAIN GREY** Stable




**CHERRY** Warm




NAVY BLUE Classic



**CAMEL** Neutral



**RED** Dynamic




**AZURE BLUE** Serenity




FOREST GREEN Natural



BLACK Mystery



#### SAKURA PINK Optimism



#### **PURPLE** Inspiring



#### AQUA GREEN Nurturing

## **14** FABRICS AND PATTERNS **VIBRANT COLORS THAT LIFT MOODS**

**UNITEX**<sup>™</sup>, is antibacterial, comfortable & stain resistant. It is Invented as a soft & anti-wrinkle covering for all radiation protection apparel.




CAMOUFLAGE



LEOPARD




BAMBINO



RETRO



BUTTERFLY



PUPPIES



HOOP





Available in various sizes Available in 0.25, 0.35 mm Pb

#### **PANORAMIC APRON**

Perfect for use with standing exposure on panoramic or tele radiographic and for digital volume tomography.



#### **DENTAL APRON**


Total protection for the patient in dentist's chair

# **16** INNER FABRIC

#### UNITEX<sup>™</sup> 360°

A futuristic range of apparel lining, made from breakthrough bamboo charcoal fiber. Organic in DNA, cutting edge in process, this new-age fabric scores high on function and fashion.

#### **THE PROCESS**



### **OUTLAST®- A UNIQUE TEMPERATURE REGULATED FABRIC**

Outlast<sup>®</sup> technology uses phase change materials (PCM) to absorb, store and release heat to achieve an optimum heat and moisture management. This gives Outlast® technology the ability to continually regulate the skin's micro climate. As the skin gets hot, the heat is absorbed, and as it cools, that heat is released. Outlast® technology does not only balance temperature fluctuations but reduces moisture.

#### THE EFFECT

- PCM changes from solid to liquid and back
- Absorbs body heat
- Stores body heat
- Releases body heat
- Maintains comfort zone

2	
Heat is stored	

2


#### 1 Excess heat is absorbed

Cooling effect

#### THE BENEFITS

- Less overheating
- Less chilling
- Less perspiration
- More comfort
- Not too hot, not too cold, just right.

Studies show that Outlast® products can reduce the sweat production average between 30% to 50%





# **17**EYE PROTECTION

### **CLEAR VISION**

European Standards 2013/59/Euratom has reduced the annual eye lens dose limit from 150 mSv/ year to **20 mSv/year** for medical staff.

Available with 0.75mm Pb frontal and 0.5mm Pb lateral protection

• Optical corrections are available as single-vision, bifocal and progressive lenses.



WITHOUT EYE PROTECTION



WITH EYE PROTECTION



### BRAZEN

bridge for comfort

#### RABID

maximum protection

OAKLEY **OAKLEY EYE WEAR** 

### GASCAN

Light-weight wraparound frame that is comfortable for all day use

Fitting frames with additional nose

Frames designed to provide clarity and



#### **FRONT & SIDE**

The user from frontal as well as lateral radiation



#### **MAXX 30**

Sporty model with padding on temples & bridge of the nose



#### **FIT OVER**

A universal solution to fit over the user's eye-wear with complete clarity and protection



#### VIENNA

Transparent design contrasted with blue, this x-ray protective glass has special material offering lateral protection

#### SYDNEY

of the eye

#### ΤΟΚΥΟ



### SPORTY

helps protect the complete eye, with no danger of lateral radiation

Metal frame in silver color with adjustable silicone nose pads, this eye-wear has an all round protection

The eye-wear has been designed to protect the eyes from scattered radiation from all angles of incidence

Wraparound style of these frames

1505



#### ROME

Rectangular style frame with high quality permanent side shields



#### PARIS

Wraparound style with a rectangular shape and a durable and light-weight design



#### BERLIN

Sporty, light-weight frame with integrated lateral protection



#### HELSINKI

Wraparound style and a light-weight design featuring removable side shields







#### LONDON

Wraparound style and a light-weight design featuring removable side shields

### **NEW YORK**

Wraparound style and a lightweight design featuring removable side shields with temple grip and comfortable nose pads

#### **UNIVERSAL SIDE SHIELDS**

Innovative lateral eye radiation protection with optimum hygiene

## **8** PRESCRIPTION SERVICE **ENHANCED VISION**

Prescription service is available in single vison, bifocal and progressive versions. A template for Prescription Form for our Radiation Protection Glasses with corrective lenses is available.



Date: \_\_\_

### **PRESCRIPTION FORM**

Prescription Values for Radiation Protection Glasses with Corrective Lenses

Customer/Name: \_\_\_\_\_ Date of Birth: \_\_\_\_\_

#### Please Tick: 🗸

Single Vision: Lined Bifocal: Progressive Bifocal: Anti-Reflective coating: / Fog freecoating: Eyewear name and/or model number:

SPHERICAL (SPH.) CYLINDRICAL (CYL.) AXIS (A) ADDITION (ADD) OD (oculus dextrus) | Right eye OD (oculus sinister) Left eye Right eye PD (pupillary distance) Left eye

Comments if any: \_\_\_\_\_

\*Please attach doctor's prescription


Signature:

# **19**CT SHIELD

Managing radiation doses during CT scans is key to avoid potential cancer risks and radiation damage to the sensitive tissues of the breast, thyroid gland and eyes.

#### **EYE SHIELD**

Reduces over 60% radiation to the eye from the CT scanner



#### **BREAST SHIELD**

Reduces cancer risk during CT scan by providing over 67% protection to the sensitive tissues of the breast

#### **THYROID SHIELD**

Provides over 65% protection to the sensitive thyroid gland and prevents tissue damage



#### **PEDIATRIC SHIELD**

Protects children who are more sensitive to radiation than adults

## **20** HAND PROTECTION **PROTECT HANDS FROM SCATTERED RADIATION**


Each disposable glove is manufactured in our bacteria-free, state of the art dipping facility.



#### ALL OUR GLOVES ARE GAMMA STERILIZED AND FEATURE:


- A lead-free proprietary hypoallergenic formulation
- Maximum contact and tactile sensitivity
- Textured surface for perfect grip
- 100% sterile and powder-free
- Elimination of finger fatigue with a unique design
- Ease of wear and removal with rolled and beaded cuff
- Latex or Latex free

RADIATION ATTENUATION - Tested as per IEC 61331-1:2014 - (BBG Method)						
	LATEX GLOVES		LATEX FREE GLOVES			
DIRECT BEAM ENERGY LEVEL	PRIME	SHEER	PRIME			
60 KV	69.6%	46.8%	63.1 %			
80 KV	60.9%	39.2%	54.4 %			
100 KV	53.8%	33.5%	47.4 %			
120 KV	48.2%	29.0%	42.0 %			
LE, mm Pb 60 to 120 kV	0.047	0.024	0.039			



## **21** SHIELDS & DRAPES **RADIATION PROTECTION AT SOURCE**

RadBlock provides protection at source when placed on the patient by creating a shade around the surgeon. Protection ranges between 75% to 96% of different types. RadBlock is available in a variety of protection levels. The BrainSaver cap provides upto 90% protection.




**RADBLOCK - FEMORAL ENTRY SHIELD 5300** 




**RADBLOCK - SUBCLAVIAN** SHIELD 5400



**RADBLOCK - PERIPHERAL** SHIELD 5110




**RADBLOCK - MULTIPURPOSE** SHIELD 5511



**RADBLOCK - BRAIN SAVER** 9100

### SHADE WITH RADBLOCK DRAPES AND COVERS



APPLICATION	PROCED
INTERVENTIONAL RADIOLOGY	Dialysis
	Liver Bio
	Chemo lı
CARDIAC CATH LAB	Biliary Dr
	Vertebro
	Balloon A
	Cardiac (
	Rotationa
ENDOVASCULAR SURGERY	Trans cat
	Carotid a
	Endovas
ELECTROPHYSIOLOGY (EP) LAB	Cardiac N
	Cardiac A
	Device In

#### DURE

opsy Immobilisation Drainage oplasty Angioplasty Catheterization nal Atherectomy

theter Aortic Valve Replacement

angioplasty and Stenting

scular Aneurysm Repair

Mapping

Ablation

mplantation

# **22**STORAGE SYSTEMS

### **COMPACT AND EFFICIENT**

Our perfectly crafted range of storage system

#### **MOBILE HANGERS**

- Smart stainless steel system with brakes
- Detachable hooks for multiple uses

#### MOUNTED WALL RACKS

- Stainless steel swivel hangers for ease of use
- Space saver racks in 1,3,5 or other types are available





**1 HANGER WALL SYSTEM** 



**5 HANGER WALL SYSTEM** 

**3 HANGER WALL SYSTEM** 




LOWER BODY METALLIC **RING STORAGE** 



# **23**OVERHEAD & LOWER BODY SHIELDS

### SAFE AND DURABLE

Designed with strong, durable and overlapping lead panels 0.5mm Pb. Multiples designs and custom shields can be made for every type of equipment. All shields can be fabricated in any length or width.



#### **OVERHEAD SHIELD**

Panels that overlap with a fold-up center part

- A fold-up center part
- 70cm (height) X 100cm (width)




#### TABLE SHIELD

Shields can be swiveled or moved laterally

 Flexible curtains 70cm (height) X 92cm (width) are used as shields to fit over C-arm equipment.

#### LOWER BODY SHIELD

- Flexible overlapping panels
- Each panel is 70cm (height) X 20cm (width)



anels



# **24** TESTING & APPROVALS

UniRay x-ray protective products conform to stringent quality norms of REGULATION (EU) 2016/425 on personal protective equipment (PPE).

The quality systems at UniRay are assessed and certified as per the guidelines of EN ISO 13485:2016 and products are tested as per strict guidelines of IEC / EN61331-1:2014.

Every UniRay apron has a four-digit registration number issued by the authorized notified body SGS-Fimko (Finland).

#### User checklist for compliance of the IEC guidelines:

- Testing done in any method other than BBG is invalid
- CE certificate (Module D) as per PPE Regulation (EU) 2016/425
- The product label should show the area weight "Ws" as mentioned in the test report.
- The product label should mention compliance with tube voltage "60 kV - 150 kV"





# JniRay

Do connect with us and get more information about our other Radiation Protection Products.

 $\mathbf{n}$ 

UniRay Medical LLP Mumbai, India E: <u>sales@uniraymedical.com</u> | W: <u>www.uniraymedical.com</u>